




# STAIR-TRAC

## PORTABLE WHEELCHAIR LIFT



Powered to carry over 250 lbs. up approximately 30 floors

The Garaventa Stair-Trac is an easy to use, portable wheelchair lift designed to attach under most manual wheelchairs. Stair-Trac's powerful motor allows a small attendant to safely transport a larger passenger up and down stairways. Stair-Trac is suitable for indoor or outdoor use in public buildings and private homes. It is an affordable and reliable portable access solution.



# STAIR-TRAC

Need Design Assistance?

Call your local Garaventa Lift representative or our Design Hotline at 1.800.663.6556

## Features


- Attaches to most standard wheelchairs with push handle (adult/child)
- Powered to carry over 250 lbs. up approximately 30 floors
- Use on virtually any indoor or outdoor stairway
- Stable and self-supporting
- Durable tracks grip stairs securely
- Detachable support handle for compact storage and easy transport
- Quick-release safety belt
- Key switch restricts use to authorized persons
- Rechargeable batteries and charger
- Portable inclined platform lift
- Suitable for indoor and outdoor use
- Accommodates all types of wheelchairs

## Benefits

- Economical compared to stairway lifts and elevators
- It is portable so it can be used on various stairways within a building
- Can be used when electrical power fails or is not available
- Easily stores in a secure area when not in use
- Requires minimal maintenance
- Available immediately
- Can be used for building evacuation

## Clearances

Wheelchair loading and unloading requires a minimum straight line clearance of approximately 3m/10ft. This diagram shows typical clearances required on turnback stairways. Dimensions vary with stairway width and operator experience. Consult Garaventa or your local representative for clearance verification.



Attach the Stair-Trac to the wheelchair, recline the passenger and move them to the stairway.



Ascend or descend stairs with the press of a button. The passenger sits in a comfortable, upright position.



Stair-Trac will climb virtually any stair tread material.

## Contact us today:

**(800) 663-6556 | garaventlift.com**

© Garaventa Lift. As we are continuously improving our products, specifications outlined in this brochure are subject to change without notice. Due to the nature of the printing process, colors shown in the brochure may not be an accurate representation of the actual product.

Authorized Garaventa Lift Representative